



# Product Data Sheet



SOFT NITRILE PREMIUM GLOVES, LATEX - FREE, POWDER - FREE

Page 1 of 1

## GENERAL INFORMATION

VIOLET-BLUE	WHITE	GREEN	PINK	PRODUCT DESCRIPTION		C = MP Class R = Rule	UMDNS	Personal protection Cat.	Notified Body	Photo
7900	7900W	7900G	7900P	XS	Medical gloves for single use, intended to protect both the wearer and the patient from cross-contamination caused by microorganisms.	CI - R1, R5	11882	III	CE	
8000	8000W	8000G	8000P	S						
8100	8100W	8100G	8100P	M						
8200	8200W	8200G	8200P	L						
8300	8300W	8300G	8300P	XL						

Physical Dimensions		Physical Properties			Physical Dimensions		Physical Properties		
Glove Length (mm)	≥ 230	Test	Before Aging	After Aging	Median glove length (mm)	≥ 240	Test	Before Aging	After Aging
Palm Thickness (mm)	0.07 ± 0.02	Tensile strength (MPa)	≥ 18	≥ 16	Median palm thickness (mm)	0.07 ± 0.02	Median Force at break (N)	≥ 6	≥ 6
Finger Thickness (mm)	0.09 ± 0.02	Elongation (%)	≥ 500	≥ 400	Median finger thickness (mm)	0.09 ± 0.02			

Regulatory Compliance:	MDR 2017/745, MDD 93/42/EEC
Standards:	EN 455 1, 2, 3, 4, EN 374 1, 2, 3, 4, 5 and EN 420 + A1
Classification:	Class I (MDR 2017/745), Class I (MDD 93/42/EEC), Category III (PPE 89/686/EEC)
Application Settings:	Low risk - medical, dental, patient care, aged care, pathology lab and food handling

### Headquarters:

Comforties.com Ltd.

Lyuben Karavelov Str. 6  
9002 Varna

**E** info@comforties.com

**W** www.comforties.com

**comforties.com**