

PRODUCT SHEET



from €
499,00

Made
in
Germany

Product	Out and indoor hand disinfectant column with foot-pump and pedestal
Capacity	5 liter standard version / childs version with 3 liter
Article Nr.	CAP-med-701 available in different colors
Packing unit (VPE)	1/1, units can be branded individually
USP	A weather-resistant, powder-coated hand disinfection column with stand (can be set up freely), a high-performance industrial foot pump technology and a lockable 5 liter disinfectant tank. The device does not require a power supply. Not suitable for hand gel disinfectants.
Production	Handmade with the highest qualification and a European Safety certificate Glaso Sp.zo.o licence Nr. WA-01062020-EVB-2140, Eco certificate GlasoSp.zo.o licence Nr. WA-01062020-EKO-156.
Protection class	99,9% bacterial filter efficiency with BIOANTICEPTICA hand disinfectant
Application	Through pressure of the foot on the foot pump, the overpressure produces the liquid to exit on the top of the unit through the valve.
Application range	Suitable for outdoor and entrance areas at: - Sport facilities, arenas, studios, clubs - Food retailing - Trade fairs and events - Public facilities as schools, departments, libraries etc. - Hospitals and nursing homes
Size / weight in kg	About 11,5 kg, height 1040 mm, pedestal 300 x 300 mm

CAPULUS assumes no guarantee / liability for topicality, correctness and completeness of the information and data provided by the manufacturer.
The product may differ from the above picture.

Contact / Orders:

Tel: +49 40 307 772 15
E-Mail: info@capulus-services.de
Request: Order Form Business Customers
Internet: CAPULUS Medical Services



ALSO AVAILABLE AS A PACKAGE DEAL WITH OUR
HAND-DISINFECTANT BIOANTISEPTICA VIRUCIDAL PLUS WITH 80% BIO-ETHANOL



CE Certified hand disinfectant made from high-purity bioethanol in accordance with EG-VO 1907/2006 (REACH) and EU-VO 4527/2010 according to the WHO formula. Optimal skin tolerance thanks to organic lipid replenishing ingredient. Manufacturer certified according to ISO 13485. Proven effectiveness, removes 99.9% bacteria, viruses and fungi.

More information at:
www.capulus-services.de/shop