

# PRODUCT SHEET



<b>product</b>	FFP2 UNICO4 CAP respirator health mask (disposable)
<b>article no.</b>	cap-med-014-2
<b>packing unit</b>	02/02 pcs package / 1.350 pcs per carton
	1.350 pcs minimum order
	goods are available immediately from our warehouse in Hamburg
<b>colour</b>	white
<b>certification</b>	DIN CE/EN 149:2001 + A1:2009 FFP
	CE Body Notified Identification Nr. 2163
<b>protection class</b>	FFP2, BFE 95% high bacterial filter efficiency
<b>range of application</b>	bacteria and viruses, airborne dusts, allergens, spores and skin scales

CAPULUS assumes no guarantee / liability for topically, correctness and completeness of the information and data provided by the manufacturer.

The product may differ from the above picture.

**Contact /Order:**

Phone: +49 40 307 772 15  
 E-Mail: info@capulus-medical.com  
 Request: Inquiry Form Business Customer  
 Internet: www.capulus-medical.com

# CERTIFICATION

UNIVERSAL



**UNIVERSAL**  
CERTIFICATION

**NB 2163**

## EU TYPE EXAMINATION CERTIFICATE

**Certificate No: 2163-PPE-1329**

Respiratory protective devices, filtering half masks to protect against particles manufactured by  
[Redacted]  
102, B [Redacted]

are tested and evaluated according to

**EN 149:2001 + A1:2009 Respiratory Protective Devices -  
Filtering Half Masks to Protect Against Particles -  
Requirements, Testing, Marking**

Based on the type examination conducted with the evaluation of test reports, technical file according to Personal Protective Equipment Regulation (EU) 2016/425 Annex 5, it is approved that the product meets the requirements of the regulation.

**Product Definition**  
Model: KN95 YM02  
Filtering half mask  
Classification: FFP2 NR

Here by the manufacturer is allowed to use notified body number (2163) and can fix CE mark, as shown below, on the Category III product models given above, with:

- Issuing an appropriate EU Declaration of Conformity according to **Personal Protective Equipment Regulation (EU) 2016/425 Annex 9**.
- Ongoing successful performance in fulfilment of the requirements set out in **Personal Protective Equipment Regulation (EU) 2016/425** and harmonised standards, ensured by assessments based on **Annex 7 (Module C2)** or **Annex 8 (Module D)** of the regulation no later than 1 year from the beginning of serial production

This certificate is initially issued on **21/08/2020** and will be valid for 5 years, if there is no change in the relevant harmonised standard affecting the essential health and safety requirements.



**2163**



**Suat KACMAZ**  
UNIVERSAL CERTIFICATION  
Director

Verify the validity with the QR code

