

PRODUCT SHEET



product	FFP2 UNICO2 CAP respirator health mask (disposable)
article no.	cap-med-s-009
packing unit	3x10 vacuum packed / 30 pcs box / 1.200 pcs. per carton
	1.200 pcs minimum order
	goods are available immediately from our warehouse in Hamburg
colour	white
certification	DIN CE/EN 149:2001 FFP NR
	CE Body Notified Identification Nr. 2163
protection class	FFP2, BFE 95% high bacterial filter efficiency
range of application	bacteria and viruses, airborne dusts, allergens, spores and skin scales

CAPULUS assumes no guarantee / liability for topically, correctness and completeness of the information and data provided by the manufacturer.


The product may differ from the above picture.

Contact /Order:

Phone: +49 40 307 772 15
 E-Mail: info@capulus-medical.de
 Request: Inquiry Form Business Customer
 Internet: www.capulus-medical.de


Certification >> click here <<

FILTERING




SMS Five-tier Structure

- ◆**Non-woven surface :**
Block larger particles;
- ◆**High-efficiency filter layer :**
Effective filtering of particles ≥ 0.075 microns
- ◆**High-efficiency filter layer :**
Effective filtering of particles ≥ 0.075 microns
- ◆**Non-woven inner layer :**
Skin friendly, comfortable and breathable
- ◆**Non-woven inner layer :**
Skin friendly, comfortable and breathable




SAMDING

3D Stereo Surround Design




Low expiratory resistance, efficient filtration
fluid resistant to splash and spatter of
blood and other infectious material



Five Layer Filtration

Blocking harmful substances layer by layer



- Non-woven
- Melt blown material
- Melt blown material
- Non-woven
- Non-woven

CERTIFICATION

UNIVERSAL



UNIVERSAL
CERTIFICATION

NB 2163

CERTIFICATE OF CONFORMANCE

Certificate Nr: 2163-PPE-638/01

Respirator manufactured by
SAMDING CRAFTWORK CO., LTD

SEAT E, SECOND INDUSTRY ZONE, JINNIU, BEIYONG, DAOJIAO TOWN,
DONGGUAN CITY, 523170, GUANGDONG, CHINA.523170

Continues to fulfil the requirements of

**EN 149:2001 + A1:2009 Respiratory Protective Devices -
Filtering Half Masks to Protect Against Particles -
Requirements, Testing, Marking**

Based on the evaluation of test reports and internal quality control audit reports according to EN 149+A1:2009 and Personal Protective Equipment Regulation (EU) 2016/425 Annex VII (Module C2). This certificate implies that the manufactured products show below are in conformance with the approved EU Type Examination model and meets the requirements of the regulation. The details of compliance is given in technical report numbered 2163-PPE-639/01

Product Definition

Model	Class	EU Type Examination Certificate		
		Serial Nr.	Date	Issuing NB Nr.
B13086	FFP2	2163-PPE-638	28/04/2020	2163

Here by the manufacturer is allowed to use notified body number (2163) and can fix CE mark, as shown below, on the Category III product models given above, with;

- Issuing an appropriate EU Declaration of Conformity according to **Personal Protective Equipment Regulation (EU) 2016/425 Annex 9.**
- Taking all measures necessary so that the manufacturing process and its monitoring ensure the homogeneity of production and conformity of the manufactured PPE with the type described in the EU type examination certificate.

This certificate is issued on **28/04/2020** and will be valid for one year, until **27/04/2021** if the manufacturer makes no major change in the product designs and manufacturing processes affecting the product performance on the essential health and safety requirement.



2163



Suat KACMAZ
UNIVERSAL CERTIFICATION
Director

The validity of this certificate can be verified online.



Necip Fazil Bulvarı Keyap Sitesi E2 Blok No:44/84 Yukarı Dudullu Ümraniye - ISTANBUL - TURKEY T:+90 216 455 80 80
UNIVERSAL CERT.COM

CERTIFICATION

UNIVERSAL



UNIVERSAL
CERTIFICATION

NB 2163

EU TYPE EXAMINATION CERTIFICATE

Certificate Nr: 2163-PPE-638

Respiratory protective devices, filtering half masks to protect against particles manufactured by

SEKİRCİLER KARAGÖZLÜKÜ YOLU, İZMİR, TÜRKİYE

are tested and evaluated according to

EN 149:2001+A1:2009 Respiratory Protective Devices - Filtering Half Masks To Protect Against Particles - Requirements, Testing, Marking

Based on the type examination conducted with the evaluation of test reports, technical file according to Personal Protective Equipment Regulation (EU) 2016/425 Annex 5, it is approved that the product meets the requirements of the regulation. The details of essential requirement compliance is given in technical report numbered 2163-PPE-639.

Product Definition

Brand Name: SAMDING **Model:** B13086

Filtering half mask

Total Inwards Leakage: Class – FFP2

Here by the manufacturer is allowed to use notified body number (2163) and can fix CE mark, as shown below, on the Category III product models given above, with;

- Issuing an appropriate EU Declaration of Conformity according to **Personal Protective Equipment Regulation (EU) 2016/425 Annex 9**.
- Ongoing successful performance in fulfilment of the requirements set out in Personal Protective Equipment Regulation (EU) 2016/425 and harmonised standards, ensured by assessments based on **Annex 7 (Module C2) or Annex 8 (Module D)** of the regulation no later than 1 year from the beginning of serial production

This certificate is initially issued on **28/04/2020** and will be valid for 5 years if there is no change in the relevant harmonised standard affecting the essential health and safety requirements.



2163



Suat KACMAZ
UNIVERSAL CERTIFICATION
Director

The validity of this certificate can be verified online.

UNIVERSALCERT.COM



Necip Fazıl Bulvarı Keyap Sitesi E2 Blok No:44/84 Yukarı Dudulu Ümraniye - İSTANBUL - TURKEY T:+90 216 455 80 80