

PRODUCT SHEET



Product	MaiMed nitrile gloves solution 100 powder free
Article No.	cap-med-902
Packing Unit	100 pcs. in box / 10 boxes in carton (packaging in German)
Thickness	Cuff min. 0,05, hand min. 0,07 mm, finger 0,09
Material	Powder free nitrile gloves free of thiurams, mercaptobenzothiazoles and zink-2 mercaptobenzothiazoles
Colour / Sizes	Blue – M / L / XL
Certification, Production	EN 374-1:2016, EN 16523-1:2015 (former 374-3:2003) EN 374-4:2013, EN 374-5:2016 A (bacteria and fungi / not viruses) EN 455 1-4, ASTM D6319, ISO 9001:2015, EN ISO 13485-2016 certified
Protection class	EN 2016/425 Category III, AQL 1,5
Application range	Food Industry, Inpatient and Outpatient Care Facilities, Primary care and Cleaning Services.

CAPULUS assumes no guarantee / liability, correctness and completeness of the information and data provided by the manufacturer. The product may differ from the above picture.

Contact / Order:

Phone: +49 40 307 772 15

Product Sheet Manufacturer >> [click here](#) <<

E-Mail: info@capulus-medical.de

Request: Inquiry Form Business Customer

Internet: www.capulus-medical.de