

## PRODUCT SHEET



|                              |   |
|------------------------------|---|
| <b>product</b>               | KN95 respirator health mask (like FFP2)   |
| <b>article no.</b>           | CAP-med-104   |
| <b>packing unit</b>          | 50/50 pcs box / 2.000 pcs per cardboard box   |
| <b>colour</b>                | white or grey   |
| <b>certification</b>         | DIN CE /EN 149:2001 + A1:2009, GB 2626 – 2006<br>CE Body Notified identification Nr. 2834 |
| <b>protection class</b>      | KN95 (like FFP2), 95% high bacterial filter efficiency                                    |
| <b>range of applications</b> | bacteria and virusses<br>airborne dusts<br>allergens<br>spores<br>skin scales             |

CAPULUS assumes no guarantee / liability for topicality, correctness and completeness of the information and data provided by the manufacturer. The product may differ from the above picture.

**Contact / Orders:**

Phone: +49 40 307 772 15

Certification / Packaging Measurements >> click here <<

E-Mail: info@capulus-services.de

Request: Inquiry Form Business Customer

Internet: CAPULUS Medical Services