

PRODUCT SHEET



product	medical face mask type IIR (disposable), non-woven 3 layer, with elastic ear loops
article no.	cap-med-108
packing unit	50/50 pcs box / 2.000 pcs per cardboard box
colour	blue
certification	CE /EN 14683:2019 + AC:2019, EU medical device regulation 93/42/EEC regarding Annex VII and Annex IX certified
protection class	98% high bacterial filter efficiency
range of applications	all applications in the surgery and ward areas

CAPULUS assumes no guarantee / liability for topicality, correctness and completeness of the information and data provided by the manufacturer. The product may differ from the above picture.

Contact / Orders:

Phone: +49 40 307 772 15

[Certification / Packaging Measurements >> click here <<](#)

E-Mail: info@capulus-medical.de

Request: Inquiry Form Business Customer

Internet: www.capulus-medical.com