


PRODUCT SHEET



Product	10 L canister, Bioantiseptica-F surface disinfectant
Container	10 L in canister
Item number	CAP-med-302
Packaging unit	10 L individual packaging
Production	According to EG-VO 1907/2006 (REACH) and EU-VO 4527/2010, manufacturer certified according to ISO 13485
Protection class	99,99% bacteria, virus and spores efficiency
Product composition	Bioantiseptica-F surface disinfectant is a hydro-alcohol solution. Ingredients: 73,0 g/100g Ethanol. (80% v/v)
Utilization	To disinfect surfaces and objects spray from a distance of 30 cm. After exposure time, wipe with disposable cloth. Alternatively, spray with Bioantiseptica-F into a disposable cloth and rub directly the surfaces and objects.
Range of application	For quick disinfection and cleaning of alcohol-resistant medical inventory and surfaces of all kinds. Also suitable for particularly inaccessible places.
Dimensions / weight in g	310 x 230 x 190 (HxWxD), 9kg

CAPULUS assumes no guarantee / liability for topicality, correctness and completeness of the information and data provided by the manufacturer. The product may differ from the above picture.



Contact / Orders:

Phone: +49 40 307 772 15

Safety data sheet >> [click here](#) <<

E-Mail: info@capulus-services.de

Instructions >> [click here](#) <<

Request: Inquiry Form Business Customer

Internet: CAPULUS Medical Services