

PRODUCT SHEET



product	infrared clinical thermometer / high standard German sensor
article no.	CAP-med-601
packing unit	Individually packed, 100 pieces in a carton
package size	185 x 46 x 46mm
weight & measurements	85g without battery / 170 x 45 x 45mm
battery	2 * AAA
certification	DIN ISO 13485 – CE Notified Body Nr. 0123 – IP22
colour	white
scope	forehead & ear
features	<ul style="list-style-type: none"> - 1 second fast measurements / 35 sets memories - LCD backlit fever alarm: yellow light from 37,4°C (from 99,3°F) red light from 38°C (from 100,4°F) - mute button - automatic shutdown after 10 seconds after measurement
operational conditions	temperature: 10°C – 40°C (50°F – 104°F) humidity: 15 -95% RH, non-condensing atmospheric pressure: 86 – 106 kPA
measuring range forehead & ear	32,0°C – 42,9°C (89,6°F – 109,2°F)
display units	0,1°C / °F

CAPULUS assumes no guarantee / liability for topicality, correctness and completeness of the information and data provided by the manufacturer. The product may differ from the above picture.

Contact / Orders:

Phone: +49 40 307 772 15

E-Mail: info@capulus-services.de

Request: Inquiry Form Business Customer

Internet: CAPULUS Medical Services

