

PRODUKTDATENBLATT



Produkt	FFP2 UNICO4 CAP Atemschutzmaske (Einweg)
Artikelnr.	cap-med-014
Verpackungseinheit (VPE)	10/10 Stk. Box / 1.350 Stk. im Karton (Verpackung in Englisch)
Beschreibung	FFP2 Einweg Atemschutzmaske ohne Filter, flach faltbar in der Mitte, mit elastischen Ohrschlaufen
Farbe	Weiß
Zertifizierung	DIN CE/EN 149:2001 + A1:2009 FFP2 CE Body Notified Identifikations-Nr. 2163
Schutzklasse	FFP2, BFE 95% hohe Bakterienfilterung
Anwendungsbereich	Bakterien und Viren, luftgetragene Stäube, Allergene, Sporen und Hautschuppen

CAPULUS übernimmt keine Gewährleistung/ Haftung für die Aktualität, Richtigkeit und Vollständigkeit der vom Hersteller bereitgestellten Informationen und Daten. Das Produkt kann vom o.g. Bild abweichen.

Kontakt / Bestellungen:

Tel: +49 40 307 772 15

Zertifizierung >> hier anklicken <<

E-Mail: info@capulus-medical.de

Anfragen: Bestellformular Geschäftskunden

Internet: www.capulus-medical.de

ZERTIFIZIERUNG

UNIVERSAL


UNIVERSAL
CERTIFICATION
NB 2163

EU TYPE EXAMINATION CERTIFICATE

Certificate No: 2163-PPE-1329

Respiratory protective devices, filtering half masks to protect against particles manufactured by
[Redacted]
102, B [Redacted]

are tested and evaluated according to

**EN 149:2001 + A1:2009 Respiratory Protective Devices -
Filtering Half Masks to Protect Against Particles -
Requirements, Testing, Marking**

Based on the type examination conducted with the evaluation of test reports, technical file according to Personal Protective Equipment Regulation (EU) 2016/425 Annex 5, it is approved that the product meets the requirements of the regulation.

Product Definition
Model: KN95 YM02
Filtering half mask
Classification: FFP2 NR

Here by the manufacturer is allowed to use notified body number (2163) and can fix CE mark, as shown below, on the Category III product models given above, with:

- Issuing an appropriate EU Declaration of Conformity according to **Personal Protective Equipment Regulation (EU) 2016/425 Annex 9**.
- Ongoing successful performance in fulfilment of the requirements set out in **Personal Protective Equipment Regulation (EU) 2016/425** and harmonised standards, ensured by assessments based on **Annex 7 (Module C2)** or **Annex 8 (Module D)** of the regulation no later than 1 year from the beginning of serial production

This certificate is initially issued on **21/08/2020** and will be valid for 5 years, if there is no change in the relevant harmonised standard affecting the essential health and safety requirements.


2163


Suat KACMAZ
UNIVERSAL CERTIFICATION
Director

Verify the validity with the QR code
