

# PRODUKTDATENBLATT



**AKTIONSPREIS**  
GÜLTIG BIS 30. OKT

<b>Produkt</b>	FFP2 NR D UNICO2 CAP Atemschutzmaske (Einweg)
<b>Artikelnr.</b>	cap-med-s-009
<b>Verpackungseinheit (VPE)</b>	30/30 Stck. Box / 1.200 Stck. Karton je Verpackungseinheit
	2.400 Stck. Mindestbestellmenge (Aktionspreis bis zum 30.09.2020)
	Ware ist sofort verfügbar ab unserem Lager in Hamburg
<b>Farbe</b>	weiß
<b>Zertifizierung</b>	DIN CE/EN 149:2001 + A1:2009, FFP2 NR D*
	CE Body Notified Identification Nr. 2163
<b>Schutzklasse</b>	FFP2, BFE 95% Bakterienfilterung
<b>Anwendungsbereich</b>	Bakterien und Viren, Luftgetragene Stäube, Allergene
	Sporen, Hautschuppen
	* D = haben diese die Dolomitstaubprüfung bestanden und sind daher zur Nutzung in einer sehr staubigen Umgebung geeignet.

ART.-NR.	PRODUKT	VK-PREISE € NT./STK.			
		bis 10.000	bis 25.000	bis 50.000	ab 50.000
Geschäftskunden	Abnahmemengen				
cap-med-109	FFP2 NR D / UNICO2 Atemschutzmaske, Filter 95%	1,25	1,15	1,05	0,95

CAPULUS übernimmt keine Gewährleistung/ Haftung für die Aktualität, Richtigkeit und Vollständigkeit der vom Hersteller bereitgestellten Informationen und Daten. Das Produkt kann vom o.g. Bild abweichen.

**Kontakt / Bestellungen:**

Tel: +49 40 307 772 15  
E-Mail: info@capulus-medical.de  
Internet: www.capulus-medical.de



## SMS Five-tier Structure


◆ **Non-woven surface :**  
Block larger particles;

◆ **High-efficiency filter layer :**  
Effective filtering of particles  $\geq 0.075$  microns

◆ **High-efficiency filter layer :**  
Effective filtering of particles  $\geq 0.075$  microns

◆ **Non-woven inner layer :**  
Skin friendly, comfortable and breathable

◆ **Non-woven inner layer :**  
Skin friendly, comfortable and breathable



## SAMDING

### 3D Stereo Surround Design




Low expiratory resistance, efficient filtration  
fluid resistant to splash and spatter of  
blood and other infectious material



### Five Layer Filtration

Blocking harmful substances layer by layer



# ZERTIFIZIERUNG

UNIVERSAL



**UNIVERSAL**  
CERTIFICATION

**NB 2163**

## CERTIFICATE OF CONFORMANCE

**Certificate Nr: 2163-PPE-638/01**

Respiratory protective devices, filtering half masks to protect against particles manufactured by  
  
 at the following manufacturing site  
 SI  /N,  
 DONDUĞAN ÇIĞI, 525/110, GENÇLİK BÖLÜMÜ, KATILIMCI YOLU, ÜSKÜDAR/İSTANBUL

Continues to fulfil the requirements of

**EN 149:2001 + A1:2009 Respiratory Protective Devices -  
Filtering Half Masks to Protect Against Particles -  
Requirements, Testing, Marking**

Based on the evaluation of test reports and internal quality control audit reports according to EN 149+A1:2009 and Personal Protective Equipment Regulation (EU) 2016/425 Annex VII (Module C2). This certificate implies that the manufactured products show below are in conformance with the approved EU Type Examination model and meets the requirements of the regulation. The details of compliance is given in technical report numbered 2163-PPE-639/01

**Product Definition**

Model	Class	EU Type Examination Certificate		
		Serial Nr.	Date	Issuing NB Nr.
B13086	FFP2	2163-PPE-638	28/04/2020	2163

Here by the manufacturer is allowed to use notified body number (2163) and can fix CE mark, as shown below, on the Category III product models given above, with;

- Issuing an appropriate EU Declaration of Conformity according to **Personal Protective Equipment Regulation (EU) 2016/425 Annex 9**.
- Taking all measures necessary so that the manufacturing process and its monitoring ensure the homogeneity of production and conformity of the manufactured PPE with the type described in the EU type examination certificate.

This certificate is issued on **28/04/2020** and will be valid for one year, until **27/04/2021** if the manufacturer makes no major change in the product designs and manufacturing processes affecting the product performance on the essential health and safety requirement.



2163



Suat KACMAZ  
UNIVERSAL CERTIFICATION  
Director

The validity of this certificate can be verified online.



Necip Fazıl Bulvarı Keyap Sitesi E2 Blok No:44/84 Yukarı Dudullu Ümraniye - İSTANBUL - TURKEY T:+90 216 455 80 80  
**UNIVERSAL CERT.COM**



