

PRODUKTDATENBLATT



AKTIONSPREIS
gültig bis 30. Nov

Produkt	FFP2 NR D UNICO2 CAP Atemschutzmaske (Einweg)
Artikelnr.	cap-med-s-009
Verpackungseinheit (VPE)	30/30 Stck. Box / 1.200 Stck. Karton je Verpackungseinheit
	1.200 Stck. Mindestbestellmenge
	Ware ist sofort verfügbar ab unserem Lager in Hamburg
Farbe	weiß
Zertifizierung	DIN CE/EN 149:2001 + A1:2009, FFP2 NR D*
	CE Body Notified Identification Nr. 2163
Schutzklasse	FFP2, BFE 95% Bakterienfilterung
Anwendungsbereich	Bakterien und Viren, Luftgetragene Stäube, Allergene
	Sporen, Hautschuppen
	* D = haben diese die Dolomitstaubprüfung bestanden und sind daher zur Nutzung in einer sehr staubigen Umgebung geeignet.

ART.-NR.	PRODUKT	VK-PREISE € NT./STK.			
		bis 10.000	bis 25.000	bis 50.000	ab 50.000
cap-med-109	FFP2 NR D / UNICO2 Atemschutzmaske, Filter 95%	1,45	1,35	1,25	1,15

CAPULUS übernimmt keine Gewährleistung/ Haftung für die Aktualität, Richtigkeit und Vollständigkeit der vom Hersteller bereitgestellten Informationen und Daten. Das Produkt kann vom o.g. Bild abweichen.

Kontakt / Bestellungen:

Tel: +49 40 307 772 15
E-Mail: info@capulus-medical.de
Internet: www.capulus-medical.de



SMS Five-tier Structure


◆ **Non-woven surface :**
Block larger particles;

◆ **High-efficiency filter layer :**
Effective filtering of particles ≥ 0.075 microns

◆ **High-efficiency filter layer :**
Effective filtering of particles ≥ 0.075 microns


◆ **Non-woven inner layer :**
Skin friendly, comfortable and breathable

◆ **Non-woven inner layer :**
Skin friendly, comfortable and breathable



SAMDING

3D Stereo Surround Design




Low expiratory resistance, efficient filtration
fluid resistant to splash and spatter of
blood and other infectious material



Five Layer Filtration

Blocking harmful
substances layer by layer



ZERTIFIZIERUNG

UNIVERSAL



UNIVERSAL
CERTIFICATION

NB 2163

CERTIFICATE OF CONFORMANCE

Certificate Nr: 2163-PPE-638/01

Respiratory protective devices, filtering half masks to protect against particles manufactured by

 at the following manufacturing site
 SI /N,
 DONDUĞAN ÇIĞIRI, 525/110, GENÇLİK BÖLGE, KATILIMCI YOLU, ÜSKÜDAR/İSTANBUL

Continues to fulfil the requirements of

**EN 149:2001 + A1:2009 Respiratory Protective Devices -
Filtering Half Masks to Protect Against Particles -
Requirements, Testing, Marking**

Based on the evaluation of test reports and internal quality control audit reports according to EN 149+A1:2009 and Personal Protective Equipment Regulation (EU) 2016/425 Annex VII (Module C2). This certificate implies that the manufactured products show below are in conformance with the approved EU Type Examination model and meets the requirements of the regulation. The details of compliance is given in technical report numbered 2163-PPE-639/01

Product Definition

Model	Class	EU Type Examination Certificate		
		Serial Nr.	Date	Issuing NB Nr.
B13086	FFP2	2163-PPE-638	28/04/2020	2163

Here by the manufacturer is allowed to use notified body number (2163) and can fix CE mark, as shown below, on the Category III product models given above, with;

- Issuing an appropriate EU Declaration of Conformity according to **Personal Protective Equipment Regulation (EU) 2016/425 Annex 9**.
- Taking all measures necessary so that the manufacturing process and its monitoring ensure the homogeneity of production and conformity of the manufactured PPE with the type described in the EU type examination certificate.

This certificate is issued on **28/04/2020** and will be valid for one year, until **27/04/2021** if the manufacturer makes no major change in the product designs and manufacturing processes affecting the product performance on the essential health and safety requirement.



2163



Suat KACMAZ
UNIVERSAL CERTIFICATION
Director

The validity of this certificate can be verified online.



Necip Fazıl Bulvarı Keyap Sitesi E2 Blok No:44/84 Yukarı Dudullu Ümraniye - İSTANBUL - TURKEY T:+90 216 455 80 80
UNIVERSAL CERT.COM

ZERTIFIZIERUNG

UNIVERSAL



UNIVERSAL
CERTIFICATION

NB 2163

EU TYPE EXAMINATION CERTIFICATE

Certificate Nr: 2163-PPE-638

Respiratory protective devices, filtering half masks to protect against particles manufactured by

 SE... ..

are tested and evaluated according to

EN 149:2001+A1:2009 Respiratory Protective Devices - Filtering Half Masks To Protect Against Particles - Requirements, Testing, Marking

Based on the type examination conducted with the evaluation of test reports, technical file according to Personal Protective Equipment Regulation (EU) 2016/425 Annex 5, it is approved that the product meets the requirements of the regulation. The details of essential requirement compliance is given in technical report numbered 2163-PPE-639.

Product Definition

Brand Name: SAMDING **Model:** B13086
 Filtering half mask
Total Inwards Leakage: Class – FFP2

Here by the manufacturer is allowed to use notified body number (2163) and can fix CE mark, as shown below, on the Category III product models given above, with;

- Issuing an appropriate EU Declaration of Conformity according to **Personal Protective Equipment Regulation (EU) 2016/425 Annex 9.**
- Ongoing successful performance in fulfilment of the requirements set out in Personal Protective Equipment Regulation (EU) 2016/425 and harmonised standards, ensured by assessments based on **Annex 7 (Module C2) or Annex 8 (Module D)** of the regulation no later than 1 year from the beginning of serial production

This certificate is initially issued on **28/04/2020** and will be valid for 5 years if there is no change in the relevant harmonised standard affecting the essential health and safety requirements.



2163



Suat KACMAZ
UNIVERSAL CERTIFICATION
Director

The validity of this certificate can be verified online.

Necip Fazıl Bulvarı Keyap Sitesi E2 Blok No:44/84 Yukarı Dudulu Ümraniye - İSTANBUL - TURKEY T:+90 216 455 80 80

UNIVERSALCERT.COM