

PRODUCT SHEET



Product	BeeSana SMS protective gown 40g, non-sterile, disposable, 3-layer
Article No.	cap-med-cap-med-410
Packing unit	10 pcs. per case, 5 cases per carton, altogether 50 gowns in carton (packaging in German)
Material	40g, Spunbond-nonwoven, Meltblown-nonwoven (SMS) material
Colour / Size	blue – one size fits all (115 x 150 cm)
Certification	Notified Body Identification-Nr. 3597, EN ISO 13485:2016 and EN 13795 certified
Protection class	93/42/EEC Medical product of class 1, virus and germ tightness by DIN EN 14126:2005, DIN EN ISO 16604, DIN EN ISO 22612, DIN EN ISO 22610
Qualities	Effective barrier against fluids and germs Opaque, low-lint, tear-resistant, breathable, high wearing comfort Fastening ties at neck and waist, cuffs at the end of the sleeves
Application range	Suitable for use on isolation wards, in elderly care, for hospital staff, patients and visitors
Measurements / Weight	39,5 x 29 x 31 cm (LxDxH), 9,8 kg carton

CAPULUS assumes no guarantee / liability, correctness and completeness of the information and data provided by the manufacturer. The product may differ from the above picture.

Contact / Order:

Phone: +49 40 307 772 15
 E-Mail: info@capulus-services.de
 Request: Inquiry Form Business Customer
 Internet: CAPULUS Medical Services

