

PRODUCT SHEET



Product	BeeSana PP protective visitor gown 23g, non-sterile, disposable
Article No.	cap-med-412
Packing unit	10 pcs. per case, 10 cases per carton, altogether 100 gowns in carton (packaging in German)
Material	23g, Polypropylene-nonwoven (PP) material
Colour / Size	green – one size fits all (110 x 140 cm)
Certification	Notified Body Identification-Nr. 3685, EN ISO 13485:2016 certified
Protection class	93/42/EEC Medical product of class 1
Qualities	Protective gown with fastening ties at neck and waist, sleeves come with sewn-in elastic rubber ties
Application range	Leightweight, breathable material, very tear-resistant while also moisture-repellent, non sterile protective gown for personnel and visitors
Measurements / Weight	39,5 x 26,5 x 31 cm (LxDxH), 4,5 kg carton

CAPULUS assumes no guarantee / liability, correctness and completeness of the information and data provided by the manufacturer. The product may differ from the above picture.

Contact / Order:

Phone: +49 40 307 772 15
 E-Mail: info@capulus-services.de
 Request: Inquiry Form Business Customer
 Internet: CAPULUS Medical Services